

Narrow Bandpass Filters 808-855nm

Code	Central Wavelength (nm)	Bandwidth (nm)	Peak Transmission (%)	Blocking Range (nm)	Optical Density	Diameter, mounted (mm)	Aperture	Delivery	Price
40112-001	808±2	10±2	> 90	200-1150	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-002	808±3	15±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-003	808±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-004	808±5	30±5	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
414-019	808nm	10nm	> 80	200 ~1200nm	≥OD4	6mm	4mm	Inquire	Inquire
414-020	808nm	10nm	> 80	200 ~1200nm	≥OD4	8mm	6mm	Inquire	Inquire
414-021	808nm	10nm	> 80	200 ~1200nm	≥OD4	10mm	6mm	Inquire	Inquire
414-022	808nm	10nm	> 80	200 ~1200nm	≥OD4	10mm	7mm	Inquire	Inquire
414-023	808nm	10nm	> 80	200 ~1200nm	≥OD4	12mm	8mm	Inquire	Inquire
414-024	808nm	10nm	> 80	200 ~1200nm	≥OD4	15mm	12mm	Inquire	Inquire
40112-005	810±2	11±2	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-006	810±3	20±3	> 90	200-1200	OD5	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-007	810±5	30±5	> 90	200-1200	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-008	813±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-009	813±5	30±5	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-010	815±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-011	815±5	30±5	> 90	200-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-012	815±5	35±5	> 90	200-1200	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-013	820±10	200±10	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-014	820±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-015	820±5	30±5	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-016	825±2	10±2	> 90	200-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-017	825±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-018	830±2	10±2	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire

40112-019	830±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-020	835±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-021	840±3	20±3	> 90	200-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-022	840±8	50±8	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-023	845±2	10±2	> 90	200-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-024	845±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-025	850±2	10±2	> 90	200-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-026	850±3	10±2	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-027	850±3	20±3	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-028	850±5	30±5	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-029	850±8	50±8	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-030	850±8	60±8	> 90	200-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-031	855±2	10±2	> 90	200-1200	OD5	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-032	855±5	30±5	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire
40112-033	855±8	50±8	> 90	400-1100	OD4	Φ12.7/Φ25.4	/	Inquire	Inquire